

Top 10 Training Tips



- 1** **VERY IMPORTANT:** Do not squeeze the brake when the person is entering or exiting the *staxi*® chair.
- 2** It is safer and easier to enter and exit the *staxi*® chair from the side. By lifting one armrest, the person can be seated. Once seated, the person pivots their body into a forward facing position. The armrest is lowered once they are comfortably in position.
- 3** Ensure the person is positioned comfortably and safely on the seat with both feet on the footrests.
- 4** The footrests can be in the lowered position at all times. This will encourage side entry and means less bending for the caregiver.
- 5** To push the *staxi*®, squeeze the brake bar and handle bar together and move.
- 6** To stop, slow the *staxi*® and then release the brake bar. The brake will automatically lock when the brake bar is released.
- 7** When stopped, the person can exit the *staxi*® by lifting the armrest the same way as on entry. The person can pivot to a side facing position on the seat placing their feet firmly on the ground and then stand.
- 8** To move through a closed doorway, it is easier to use one hand to pull the *staxi*® backwards through the door and the other hand to hold open the door.
- 9** *staxi*® has two convenient storage areas; one area under the seat for suitcases and larger bags as well as a second area behind the seat for magazines and documents. This second area also has hooks to hang items such as purses or smaller parcels.
- 10** The maximum weight capacity for Medical and Commercial *staxi*® is 600lbs, for Ranger Staxi 800lbs and for Max Staxi 1000lbs. The weight capacity for your chair is shown on the handle bar label. The armrests may remain in the raised position if required when transporting larger persons.



Toll Free Phone +1-877-677-8294

Email: info@staxi.com Website: www.staxi.com



Strongest Design



Hardest to Steal



Safest to Use

Patent No.(US) 6,315,306 Patent No.(US) 6,161,848 Patent No.(US) 6,412,795 Patent No.(CAN) 2250184 Patent No.(EUR) P11385GBEP Patent No.(UK) EP0889709

staxi®

Rolling out Excellence