

Recommended Pre-use Inspection Checks

Every chair is thoroughly inspected prior to leaving our facilities in accordance with STAXI[®] Corporation Ltd's stringent quality control procedures to ensure that each **staxi**[®] is in perfect working order. However, we strongly recommend that your new **staxi**[®] should now be inspected prior to use as on some occasions shipping can cause some movement to the mechanical parts of the chair. Prior to putting your **staxi**[®] into service we recommend that you undertake the following checks:

- 1 See our 'Damaged in Transit Policy' for the correct procedures to follow if any damage to either the chairs or the packing materials is visible, prior to accepting this shipment. [Click here](#) to view this in the 'Wow-Care™' section at www.staxi.com**
- 2 Visually inspect each chair to see if there are any signs of shipping damage that was not apparent prior to unpacking your new **staxi**[®] chairs.**
- 3 Remove all of the outer packaging, including the protective foam around the push bar.**
- 4 ENSURE THAT THE BRAKES ARE WORKING CORRECTLY. Follow the checks on this Pre-use Inspection Checklist and as recommended in our 'Maintenance and Inspection Schedule'. [Click here](#) to download this from the 'Quick Start' section at www.staxi.com**
- 5 Inspect the rear wheels and front casters for any signs of damage.**
- 6 Ensure that the armrests and footrests are working correctly. Follow the checks as recommended in our Maintenance and Inspection Schedule.**
- 7 When each **staxi**[®] chair passes this inspection checklist, remove the red tag and enjoy your new **staxi**[®].**

NOTE: The Inspection Checklist and further Use and Care manuals and videos are available from our website. [Click here](#) to download from the 'Quick Start' section at www.staxi.com

We strongly recommend you read and download these manuals and also watch the video presentation before your **staxi**[®] is put into service.

IMPORTANT: If you notice any damage to your **staxi**[®] following use of this checklist, please inform Wow-Care™ immediately on:

Toll Free Phone: 1-866-969-1463 Ext 207 or **Toll Free Phone:** 1-866-WOW-1463 Ext 207
or **Email:** Wow-Care@staxi.com



Strongest Design




Hardest to Steal




Safest to Use


IMPORTANT: Brake Inspection Checks



- 1 Push brake lever and roll chair forward
- 2 If brakes 'click' with lever depressed then brakes need adjusting.
- 3 Squeeze the brake bar against the push bar and hold. Adjust the brake by turning the nut at the top of black lower brake rod. While doing this watch the pin that engages in the brake locking cogs. When this pin is just clear of the cogs the brake is adjusted correctly. Release the brake bar and make certain that when released, the pin is seated in the bottom of the brake locking cogs.



- 4 If brake is not falling into socket, loosen bolts very slightly or until the brake falls. See diagram 3 (A & B)



Alternatively visit the 'Customer Info' section of our website at www.staxi.com to view the 'Monthly Maintenance' video presentation.