





| | |
|---------|-------------------|
| Product | Aisle Chair AC020 |
|---------|-------------------|

| | |
|-------------------------------------|---|
| Manufacturer | Staxi Corp. 424 Williamson- Zebulon Rd. Williamson GA 30292 USA |
| EU Authorized Representative | MedEnvoy Global B.V. Prinses Margrietplantsoen 33 – Suite 123. 2595 AM The Hague. The Netherlands |
| UK Authorized Representative | MedEnvoy UK Limited. 85 Great Portland St, First Floor London. W1W 7LT United Kingdom |

| 1.0 Classification and Intended Use | |
|---|---|
| <p>Staxi Aisle Chair is a class I medical device under Regulation (EU) 2017/745 (for MDR). EMDN Code: V080502.</p> <p>The Staxi Aisle Chair is intended to be used in the transport of non-ambulatory passengers from their personal mobility assistance equipment to their seat on an aircraft, and vice-versa, under the direction of 2 attendants. Typical use involves a single transport episode of transient duration (less than 1 hour).</p> |  <p>Fig 1. Staxi Aisle Chair</p> |

| 2.0 Chair Part Codes – incl. Colors | | | |
|-------------------------------------|---|--|--|
| AC020-X |  | | |

| 3.0 Instructions for safe use: |
|---|
| <ul style="list-style-type: none"> • To be used transporting passengers from aircraft door to seat only. • Before use, visually inspect the chair to confirm that all components are securely fastened, undamaged, and 4 wheels roll smoothly. • Ensure that the brake is functioning properly by testing that the chair will not roll forward while the brake is engaged. • Ensure the secondary safety brake bar is engaging properly. • If installed, test that the armrest does not free-fall from the raised position. Ensure all safety straps and buckles are secure and operational . • If any issues or potential issues are detected during inspection, do not use. Remove the chair from service and contact maintenance. • Safe and unobstructed access is provided via side entry, with the armrest in the raised position if installed. • Use of side entry and allows footrests to remain lowered at all times. • If transfer process calls for crossing a ledge or gap into the aircraft, use the front and rear lift bars to cross the gap. • Always use two attendants to perform transfers, one at the front and one at the rear of the transport chair. |



- Attendants are to ensure the chair does not move during the transfer
- Transfer passenger as per the responsible companies 'Transfer SOP'
- **Harness, both thigh straps, and lower leg strap must be secured before any movement.**
Do not use if any of the above is missing or defective
- To move chair:
 - Before disengaging the brake or moving the chair, the passengers' feet must be placed firmly on the lowered footrests.
 - The lower leg restraint should be firmly adjusted to ensure feet and legs stay secure.
 - Lift the safety bar then disengage the brake by squeezing the brake bar. The chair is now prepared for movement.
 - Engage the brake by releasing the brake bar. The chair must be stationary before the brake is engaged. The safety bar will engage
 - Do not disengage the brake during passenger entry or exit.
- All points on the safety harness must be used and adjusted specifically for each passenger. The thigh strap and lower leg strap must be used and adjusted to fit each passenger. After onboard passenger transfer, harness and all straps should be secured
- Weight Capacity: 500lbs / 227kg.
- Keep the chair away from open flames and sources of intense heat.
- Attendants must stay with chair when occupied.
- Do not overload the chair.
- Avoid tip-overs and falls.
 - Only one passenger at a time.
 - Do not attempt a U-turn or 360 maneuver on a slope, ramp or incline.
 - Do not stand on footrests.
 - Never wheel this chair on escalators or stairs.
- Use caution when maneuvering the chair, to avoid pinching arms or hands between armrests and walls or door edges when navigating narrow doorways or tight gaps.
 - Avoid tight spaces.
- Exclusively for use when passengers require assistance from their personal wheelchair to the aircraft seat. Designed to ease egress down the airplane aisle. Do not use for general transport from point-to-point in the airport.
- Never place children in car seats on this chair.

3.1 Accessory Instruction Information

The Following information is pertaining to Staxi Accessories fitted to the Chair.

- **Head Rest**
 - If the headrest is already in place on the Aisle Chair with a passenger to be transferred, adjust the headrest to support the passengers head. It is important to communicate with the passenger to ensure comfort and safety.
 - If the headrest is not already positioned in place on the chair the headrest can be located on the chair's lower cross-bar. Lift the T-Bar and remove the headrest. Slide the headrest mount bar into the mounting bracket and adjust for the customers comfort and safety.

4.0 Handling and Storage:

- Do not stack chairs on top of each other.
- Keep chairs indoors where possible, protected from weather to maximize lifespan.

5.0 Cleaning:

- Clean between users or when soiled using a neutral-pH, non-corrosive detergent or approved low-level disinfectant; follow the cleaning product usage label for dilution/contact time.
- Wipe only; do not immerse. Avoid solvents, abrasives, steam/pressure washing, and unapproved corrosive agents.
- After cleaning, inspect critical parts; remove from service if damaged or not functioning.

6.0 End-of-Life & Disposal Information:

The expected lifetime of the device is 10 years. This service life is based on assessment of normal wear-and-tear under typical conditions of use and routine maintenance.

Dispose of the chair and its components in accordance with local regulations. End-of-life components may be managed through normal solid-waste or recycling streams. Cushion and textile parts may be recycled where programs exist.

7.0 Declaration:**For Europe and United Kingdom**

Staxi Corp declares that the above-mentioned products meet the provisions of the following EU legislation:

- Medical Devices Regulation (EU) 2017/745
- EN ISO 15223-1:2021 for symbols
- EN ISO 20417:2021 for labelling

For the United States

This device is manufactured in accordance with the applicable requirements of the U.S. Federal Food, Drug, and Cosmetic Act and 21 CFR Part 820 Quality System Regulation.

The manufacturer is registered with the U.S. Food and Drug Administration and the device is listed as required under 21 CFR Part 807.

8.0 Notice:

If you observe or have a serious incident with the Staxi Aisle Chair, report this to Staxi and the and the competent authority of the Member State in which you are based.